



PHKG2 mouse mAb

Catalog No	YP-mAb-10245
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PHKG2
Immunogen	Synthesized peptide derived from human PHKG2 AA range: 166-216
Specificity	This antibody detects endogenous levels of human PHKG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.7.11.19;GSD9C;PHK gamma T;PHK-gamma-LT;PHK-gamma-T;PHKG2; Phosphorylase b kinase gamma catalytic chain;testis/liver isoform (PHK-gamma-T) (EC 2.7.11.19) (PSK-C3) (Phosphorylase kinase subunit gamma-2) (Serine/threonine-protein kinase PHKG2) (EC 2.7.11.1) (E; testis/liver isoform (PHK-gamma-T) (EC 2.7.11.19) (PSK-C3) (Phosphorylase kinase subunit gamma-2) (Serine/threonine-protein kinase PHKG2) (EC 2.7.11.1) (E
Observed Band	46 kDa
Calculated Molecular Weight	46 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images