



G2E3 mouse mAb

Catalog No	YP-mAb-10246
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	G2E3 KIAA1333
Immunogen	Synthesized peptide derived from human G2E3 AA range: 196-246
Specificity	This antibody detects endogenous levels of human G2E3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G2E3;G2E3 KIAA1333;G2/M phase-specific E3 ubiquitin-protein ligase (EC 6.3.2.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus . Cytoplasm . Shuttles between the nucleus and the cytoplasm. In the nucleus, delocalizes from the nucleolus to the nucleoplasm in response to DNA damage. .
Tissue Specificity	Predominantly expressed in brain, liver, kidney, testes and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images