



C144A mouse mAb

Catalog No	YP-mAb-10251
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC144A KIAA0565
Immunogen	Synthesized peptide derived from human C144A AA range: 1116-1166
Specificity	This antibody detects endogenous levels of human C144A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C144A;CCDC144A KIAA0565;Coiled-coil domain-containing protein 144A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images