



BEST3 mouse mAb

Catalog No	YP-mAb-10253
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BEST3 VMD2L3
Immunogen	Synthesized peptide derived from human BEST3 AA range: 79-129
Specificity	This antibody detects endogenous levels of human BEST3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BEST3;Bestrophin 3;Bestrophin-3;VMD2L3;BEST3 VMD2L3;Bestrophin-3 (Vitelliform macular dystrophy 2-like protein 3)
Observed Band	76-78 kDa, 50-55 kDa
Calculated Molecular Weight	668 aa, 76 kDa
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Present in skeletal muscle and weakly in brain, spinal cord, bone marrow and retina.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images