



## GRD2I mouse mAb

|                             |                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-mAb-10254                                                                                                                      |
| Isotype                     | IgG                                                                                                                               |
| Reactivity                  | Human;Mouse;Rat                                                                                                                   |
| Applications                | WB                                                                                                                                |
| Gene Name                   | GRID2IP                                                                                                                           |
| Immunogen                   | Synthesized peptide derived from human GRD2I AA range: 43-93                                                                      |
| Specificity                 | This antibody detects endogenous levels of human GRD2I                                                                            |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| Source                      | Monoclonal, Mouse,IgG                                                                                                             |
| Purification                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| Dilution                    | WB 1:500-1:2000                                                                                                                   |
| Concentration               | 1 mg/ml                                                                                                                           |
| Purity                      | ≥90%                                                                                                                              |
| Storage Stability           | -20°C/1 year                                                                                                                      |
| Synonyms                    | GRD2I;GRID2IP;Delphinin (Glutamate receptor;ionotropic;delta 2-interacting protein 1); ionotropic; delta 2-interacting protein 1) |
| Observed Band               |                                                                                                                                   |
| Calculated Molecular Weight |                                                                                                                                   |
| Cell Pathway                | Cell junction, synapse, postsynaptic cell membrane .                                                                              |
| Tissue Specificity          |                                                                                                                                   |
| Function                    |                                                                                                                                   |
| Background                  |                                                                                                                                   |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                              |
| Usage suggestions           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images