



PRS6B mouse mAb

Catalog No	YP-mAb-10255
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSMC4 MIP224 TBP7
Immunogen	Synthesized peptide derived from human PRS6B AA range: 329-379
Specificity	This antibody detects endogenous levels of human PRS6B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRS6B;PSMC4 MIP224 TBP7;26S protease regulatory subunit 6B (26S proteasome AAA-ATPase subunit RPT3) (MB67-interacting protein) (MIP224) (Proteasome 26S subunit ATPase 4) (Tat-binding protein 7) (TBP-7)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images