



ATG2B mouse mAb

Catalog No	YP-mAb-10256
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ATG2B C14orf103
Immunogen	Synthesized peptide derived from human ATG2B AA range: 266-316
Specificity	This antibody detects endogenous levels of human ATG2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf103;Autophagy-related protein 2 homolog B;ATG2B;ATG2B C14orf103
Observed Band	240-250 kDa
Calculated Molecular Weight	2078 aa, 233 kDa
Cell Pathway	Preautophagosomal structure membrane ; Peripheral membrane protein . Lipid droplet . Endoplasmic reticulum membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images