



CLD15 mouse mAb

Catalog No	YP-mAb-10261
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLDN15
Immunogen	Synthesized peptide derived from human CLD15 AA range: 66-116
Specificity	This antibody detects endogenous levels of human CLD15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLD15;CLDN15;Claudin-15
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, tight junction. Cell membrane; Multi-pass membrane protein. Tight junctions form continuous circumferential cell-cell contacts at the borders of apical and lateral cell membranes that seal the intercellular space and show up as strand-like structures in electron microscopy.
Tissue Specificity	Detected in colon (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images