



CK049 mouse mAb

Catalog No	YP-mAb-10262
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C11orf49
Immunogen	Synthesized peptide derived from human CK049 AA range: 165-215
Specificity	This antibody detects endogenous levels of human CK049
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CK049;C11orf49;UPF0705 protein C11orf49
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite . Cytoplasm, cytoskeleton . Associated with microtubules. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images