



SPTCS mouse mAb

Catalog No	YP-mAb-10263
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPG11 KIAA1840
Immunogen	Synthesized peptide derived from human SPTCS AA range: 618-668
Specificity	This antibody detects endogenous levels of human SPTCS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPTCS;SPG11 KIAA1840;Spatacsin (Colorectal carcinoma-associated protein) (Spastic paraplegia 11 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Nucleus . Cell projection, axon . Cell projection, dendrite . Mainly cytoplasmic. .
Tissue Specificity	Expressed in all structures of brain, with a high expression in cerebellum. Expressed in cortical projection neurons.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images