



PPC1A mouse mAb

Catalog No	YP-mAb-10264
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPAPDC1A DPPL2 PPAPDC1
Immunogen	Synthesized peptide derived from human PPC1A AA range: 165-215
Specificity	This antibody detects endogenous levels of human PPC1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PPC1A;PPAPDC1A DPPL2 PPAPDC1;Phosphatidate phosphatase PPAPDC1A (EC 3.1.3.4) (Phosphatidic acid phosphatase type 2 domain-containing protein 1A)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed mainly to the brain, kidney and testis, and to a lesser extent the bone marrow, thymus, prostate, liver and uterus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images