



ENTP3 mouse mAb

Catalog No	YP-mAb-10265
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ENTPD3 CD39L3
Immunogen	Synthesized peptide derived from human ENTP3 AA range: 298-348
Specificity	This antibody detects endogenous levels of human ENTP3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENTP3;ENTPD3 CD39L3;Ectonucleoside triphosphate diphosphohydrolase 3 (NTPDase 3) (EC 3.6.1.5) (CD39 antigen-like 3) (Ecto-ATP diphosphohydrolase 3) (Ecto-ATPDase 3) (Ecto-ATPase 3) (Ecto-apyrase 3) (HB6)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in adult brain, pancreas, spleen and prostate (PubMed:9676430). Moderate or low expression is seen in most tissues (PubMed:9676430). Not expressed in liver and peripheral blood leukocytes (PubMed:9676430).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images