



ULK2 Monoclonal Antibody

Catalog No	YP-mAb-10593
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ULK2 KIAA0623
Immunogen	Synthetic Peptide of ULK2
Specificity	ULK2 Monoclonal Antibody detects endogenous levels of ULK2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ULK2 KIAA0623;ULK2;unc-51-like kinase 2 (C. elegans)
Observed Band	125
Calculated Molecular Weight	112694
Cell Pathway	Cytoplasmic vesicle membrane ; Peripheral membrane protein . Localizes to pre-autophagosomal membrane.
Tissue Specificity	Brain,Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images