



Glypican-2 Monoclonal Antibody

Catalog No	YP-mAb-10865
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPC2
Immunogen	Synthesized peptide derived from human Glypican-2
Specificity	This antibody detects endogenous levels of human Glypican-2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPC2;Secreted glypican-2;FLJ38962;Glypican-2;Glypican-2 [Cleaved into: Secreted glypican-2]
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor ; Extracellular side .; [Secreted glypican-2]: Secreted, extracellular space .
Tissue Specificity	Brain,
Function	function:Cell surface proteoglycan that bears heparan sulfate. May fulfill a function related to the motile behaviors of developing neurons.,similarity:Belongs to the glypican family.,
Background	function:Cell surface proteoglycan that bears heparan sulfate. May fulfill a function related to the motile behaviors of developing neurons.,similarity:Belongs to the glypican family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images