



S38A7 mouse mAb

Catalog No	YP-mAb-11422
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC38A7 SNAT7
Immunogen	Synthesized peptide derived from human S38A7 AA range: 209-259
Specificity	This antibody detects endogenous levels of S38A7 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SNAT7;Solute carrier family 38 member 7;Sodium-coupled neutral amino acid transporter 7;S38A7;SLC38A7 SNAT7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . In neurons, located in soma and axons. .
Tissue Specificity	
Function	function:Putative sodium-dependent amino acid/proton antiporter., similarity:Belongs to the amino acid/polyamine transporter 2 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images