



YB039 mouse mAb

Catalog No	YP-mAb-11936
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LINC00116 NCRNA00116
Immunogen	Synthesized peptide derived from human YB039
Specificity	This antibody detects endogenous levels of YB039 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LINC00116 NCRNA00116;YB039
Observed Band	36850DA
Calculated Molecular Weight	15180
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein . Preferentially binds to cardiolipin relative to other common cell membrane lipids. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images