




JOS2 mouse mAb

Catalog No	YP-mAb-12248
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	JOSD2 SBBI54
Immunogen	Synthesized peptide derived from human JOS2 AA range: 62-112
Specificity	This antibody detects endogenous levels of JOS2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Cell Pathway	Cytoplasm, cytosol .
Calculated Molecular Weight	21kD
Tissue Specificity	
Function	similarity:Contains 1 Josephin domain.,
Background	This gene encodes a protein containing a Josephin domain. Josephin domain-containing proteins are deubiquitinating enzymes which catalyze the hydrolysis of the bond between the C-terminal glycine of the ubiquitin peptide and protein substrates. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4°over night