



GROβ mouse mAb

Catalog No	YP-mAb-12432
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CXCL2 GRO2 GROB MIP2A SCYB2
Immunogen	Synthesized peptide derived from human GROβ
Specificity	This antibody detects endogenous levels of Human GROβ
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CXCL2 GRO2 GROB MIP2A SCYB2;GROβ;C-X-C motif chemokine 2 (Growth-regulated protein beta;Gro-beta;Macrophage inflammatory protein 2-alpha;MIP2-alpha) [Cleaved into: GRO-beta(5-73;GRO-beta-T;Hematopoietic synergistic factor;HSF;SB-251353)]
Observed Band	21780DA
Calculated Molecular Weight	11770
Cell Pathway	Secreted.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images