



2019-covid N protein mouse mAb

Catalog No	YP-mAb-12453
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	2019-covid N protein
Immunogen	Synthesized peptide derived from 2019-covid N protein
Specificity	This antibody detects endogenous levels of 2019-covid N protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	2019-covid N protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images