



IL17F mouse mAb

Catalog No	YP-mAb-12458
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IL17F
Immunogen	Synthesized peptide derived from mouse IL17F
Specificity	This antibody detects endogenous levels of Mouse IL17F
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Il17f;interleukin 17F;IL 17F;C87042;Interleukin-17F (IL-17F;Cytokine ML-1; Interleukin-24;IL-24)
Calculated Molecular Weight	17KD
Observed Band	17-25 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

