



BDNF mouse mAb

Catalog No	YP-mAb-12464
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BDNF
Immunogen	Synthesized peptide derived from mouse BDNF
Specificity	This antibody detects endogenous levels of Mouse,Rat BDNF
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BDNF precursor form;brain derived neurotrophic factor;Brain-derived neurotrophic factor;BDNF
Observed Band 28kD	28-30 kDa
Calculated Molecular Weight	247 aa, 28 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images