



FGF7 mouse mAb

Catalog No	YP-mAb-12476
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FGF7 KGF
Immunogen	Synthesized peptide derived from human FGF7
Specificity	This antibody detects endogenous levels of Human,Rat FGF7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FGF7;FGF-7;FGF 7;Fibroblast growth factor 7;HBGF-7;FGF7 KGF;Fibroblast growth factor 7 (FGF-7;Heparin-binding growth factor 7;Keratinocyte growth factor)
Observed Band	23 kDa
Calculated Molecular Weight	22-28 kDa
Cell Pathway	Secreted.
Tissue Specificity	Epithelial cell.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images