



UPAR mouse mAb

Catalog No	YP-mAb-12481
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLAUR MO3 UPAR
Immunogen	Synthesized peptide derived from human UPAR
Specificity	This antibody detects endogenous levels of Human UPAR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLAUR MO3 UPAR;UPAR;Urokinase plasminogen activator surface receptor (UPAR;Monocyte activation antigen Mo3;CD antigen CD87)
Observed Band	118470DA
Calculated Molecular Weight	36850
Cell Pathway	Cell membrane . Cell projection, invadopodium membrane . Colocalized with FAP (seprase) preferentially at the cell surface of invadopodia membrane in a cytoskeleton-, integrin- and vitronectin-dependent manner. .; [Isoform 1]: Cell membrane ; Lipid-anchor, GPI-anchor .; [Isoform 2]: Secreted .
Tissue Specificity	Expressed in neurons of the rolandic area of the brain (at protein level). Expressed in the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images