



Podoplanin mouse mAb

Catalog No	YP-mAb-12510
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	Podoplanin
Immunogen	Synthesized peptide derived from human Podoplanin
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat Podoplanin
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T1A;GP36;GP40;Gp38;OTS8;T1A2;T11A;D2-40;T1A-2;AGGRUS;HT1A-1;PA2.26;Podoplanin
Observed Band	36 kDa
Calculated Molecular Weight	17 kDa/18 kDa/24 kDa/12 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images