



TAX1BP3 mouse mAb

Catalog No	YP-mAb-12538
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TAX1BP3 TIP1
Immunogen	Synthesized peptide derived from human TAX1BP3
Specificity	This antibody detects endogenous levels of TAX1BP3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.20% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tax interaction protein 1;Tax1 binding protein 3;TAX1BP3;TIP 1;TIP1;TAX1BP3 TIP1;Tax1-binding protein 3 (Glutaminase-interacting protein 3) (Tax interaction protein 1) (TIP-1) (Tax-interacting protein 1)
Observed Band	15-17 kDa
Calculated Molecular Weight	124 aa, 14 kDa
Cell Pathway	Cytoplasm. Nucleus. Cell membrane; Peripheral membrane protein; Cytoplasmic side. Recruited to the cell membrane by interaction with membrane proteins.
Tissue Specificity	Ubiquitous. Detected in brain, heart, kidney, lung, small intestine and skeletal muscle. Detected in various cell lines including HeLa. Weakly expressed in peripheral blood leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images