



FSTL1 mouse mAb

Catalog No	YP-mAb-12543
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FSTL1 FRP
Immunogen	Synthesized peptide derived from human FSTL1
Specificity	This antibody detects endogenous levels of FSTL1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.29% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	folliculin like 1;Folliculin like protein 1;Folliculin related protein 1;Folliculin-like protein 1;FSTL1;FSTL1 FRP;Folliculin-related protein 1 (Folliculin-like protein 1)
Observed Band	35-40 kDa
Calculated Molecular Weight	308 aa, 35 kDa
Cell Pathway	Secreted .
Tissue Specificity	Overexpressed in synovial tissues from rheumatoid arthritis (PubMed:15638044).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images