



USP28 mouse mAb

Catalog No	YP-mAb-12553
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP28 KIAA1515
Immunogen	Synthesized peptide derived from human USP28
Specificity	This antibody detects endogenous levels of USP28 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.41% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USP28 KIAA1515;USP28;Ubiquitin carboxyl-terminal hydrolase 28 (EC 3.4.19.12) (Deubiquitinating enzyme 28) (Ubiquitin thioesterase 28) (Ubiquitin-specific-processing protease 28)
Observed Band	62370DA
Calculated Molecular Weight	118470
Cell Pathway	Nucleus, nucleoplasm .
Tissue Specificity	Brain, Eye,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images