



RPL41 mouse mAb

Catalog No	YP-mAb-12562
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL41
Immunogen	Synthesized peptide derived from human RPL41 AA range: 1-25
Specificity	This antibody detects endogenous levels of RPL41 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.53% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	6S ribosomal protein L41;L41;RPL41;60S ribosomal protein L41 (HG12)
Observed Band	
Calculated Molecular Weight	3kDa; 3kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images