



HIVEP2 mouse mAb

Catalog No	YP-mAb-12563
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HIVEP2
Immunogen	Synthesized peptide derived from human HIVEP2
Specificity	This antibody detects endogenous levels of HIVEP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.54% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HIV EP2;HIV-EP2;HIVEP2;human immunodeficiency virus type I enhancer binding protein 2;Human immunodeficiency virus type I enhancer-binding protein 2;MBP 2;MBP-2;MHC binding protein 2;MHC-binding protein 2;MIBP1;SHN2; Transcription factor HIVEP2;ZAS2;ZEP2_HUMAN;ZNF4B;Transcription factor HIVEP2 (Human immunodeficiency virus type I enhancer-binding protein 2) (HIV-EP2) (MHC-binding protein 2) (MBP-2)
Observed Band	
Calculated Molecular Weight	35,39 kDa; 269kD(Calculated).
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in brain and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images