



# ISK7 mouse mAb

|                                    |                                                                                                                                                           |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-12564                                                                                                                                              |
| <b>Isotype</b>                     | IgG                                                                                                                                                       |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                           |
| <b>Applications</b>                | WB                                                                                                                                                        |
| <b>Gene Name</b>                   | SPINK7 ECG2 UNQ745/PRO1474                                                                                                                                |
| <b>Immunogen</b>                   | Synthesized peptide derived from human ISK7 AA range: 1-50                                                                                                |
| <b>Specificity</b>                 | This antibody detects endogenous levels of ISK7 at Human, Mouse,Rat                                                                                       |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, and 0.55% sodium azide.                                                                                            |
| <b>Source</b>                      | Monoclonal, mouse,IgG                                                                                                                                     |
| <b>Purification</b>                | The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.                                                  |
| <b>Dilution</b>                    | WB 1:500-2000                                                                                                                                             |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                   |
| <b>Purity</b>                      | ≥90%                                                                                                                                                      |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                              |
| <b>Synonyms</b>                    | ECG2;ECRG 2;ECRG2;ECRG-2;SPINK7;ISK7;SPINK7 ECG2 UNQ745/PRO1474;Serine protease inhibitor Kazal-type 7 (Esophagus cancer-related gene 2 protein) (ECRG-2) |
| <b>Observed Band</b>               | 20kD                                                                                                                                                      |
| <b>Calculated Molecular Weight</b> | 20kD                                                                                                                                                      |
| <b>Cell Pathway</b>                | Secreted .                                                                                                                                                |
| <b>Tissue Specificity</b>          |                                                                                                                                                           |
| <b>Function</b>                    |                                                                                                                                                           |
| <b>Background</b>                  |                                                                                                                                                           |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                      |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                         |



## Products Images

