



ZFHX4 mouse mAb

Catalog No	YP-mAb-12566
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZFHX4
Immunogen	Synthesized peptide derived from human ZFHX4
Specificity	This antibody detects endogenous levels of ZFHX4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.57% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZFH 4;ZFH4;ZFH-4;ZHF4;ZFHX4;Zinc finger homeobox protein 4 (Zinc finger homeodomain protein 4) (ZFH-4)
Observed Band 392kD	420 kDa
Calculated Molecular Weight	394 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in brain, skeletal muscle and liver. Very low expression in stomach.
Function	disease:A chromosomal aberration involving [ZFHX4] is found in one patient with ptosis. Translocation t(1;8)(p34.3;q21.12).,function:May play a role in neural and muscle differentiation (By similarity). May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 20 C2H2-type zinc fingers.,similarity:Contains 4 homeobox DNA-binding domains.,tissue specificity:Expressed in brain, skeletal muscle and liver. Very low expression in stomach.,
Background	



**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images