



Salusin- β mouse mAb

Catalog No	YP-mAb-12571
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TOR2A HEMBA1005096 PSEC0218
Immunogen	Synthesized peptide derived from human Salusin- β
Specificity	This antibody detects endogenous levels of Salusin- β at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.62% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	TORP1;Salusin-beta;Salusin-alpha;PSEC0218;Prosalusin;Salusin- β ;TOR2A HEMBA1005096 PSEC0218;Prosalusin (Torsin family 2 member A) (Torsin-2A) [Cleaved into: Salusin-alpha;Salusin-beta]
Observed Band	40-45 kDa
Calculated Molecular Weight	36 kDa
Cell Pathway	Secreted .
Tissue Specificity	Isoform 4 is ubiquitously expressed, with high level in vascular endothelial cells and vascular smooth muscle cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images