



TMEM147 mouse mAb

Catalog No	YP-mAb-12576
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM147
Immunogen	Synthesized peptide derived from human TMEM147
Specificity	This antibody detects endogenous levels of TMEM147 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.68% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NEDFLPH;NIFIE14;TMEM147;Transmembrane protein 147 (Protein NIFIE 14)
Observed Band	25kDa
Calculated Molecular Weight	25kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images