



ITPR3 mouse mAb

Catalog No	YP-mAb-12577
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ITPR3
Immunogen	Synthesized peptide derived from human ITPR3
Specificity	This antibody detects endogenous levels of ITPR3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.69% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	InsP3R3;IP3R;IP3R 3;IP3R3;ITPR3;Inositol 1;4;5-trisphosphate receptor type 3 (IP3 receptor isoform 3) (IP3R 3) (InsP3R3) (Type 3 inositol 1;5-trisphosphate receptor) (Type 3 InsP3 receptor)
Observed Band	250-300 kDa
Calculated Molecular Weight	304 kDa
Cell Pathway	Endoplasmic reticulum membrane; Multi-pass membrane protein.
Tissue Specificity	Expressed in intestinal crypt and villus epithelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images