



NDUFAF6 mouse mAb

Catalog No	YP-mAb-12581
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NDUFAF6 C8orf38
Immunogen	Synthesized peptide derived from human NDUFAF6
Specificity	This antibody detects endogenous levels of NDUFAF6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.73% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NDUFAF6;NDUFAF6 C8orf38;NADH dehydrogenase (ubiquinone) complex I; assembly factor 6 (Putative phytoene synthase); assembly factor 6 (Putative phytoene synthase)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Mitochondrion inner membrane. Peripherally localized on the matrix face of the mitochondrial inner membrane.; [Isoform 2]: Cytoplasm. Nucleus.
Tissue Specificity	Widely expressed. A lower expression is observed in lung and kidney compared to heart, muscle and liver (PubMed:27466185). In the kidney, expression is high in the basal zone of the proximal tubular cells (PubMed:27466185).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images