



LUR1 mouse mAb

Catalog No	YP-mAb-12583
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C12orf49
Immunogen	Synthesized peptide derived from human LUR1
Specificity	This antibody detects endogenous levels of LUR1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.75% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LUR1;C12orf49;UPF0454 protein C12orf49
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images