



glrx5 mouse mAb

Catalog No	YP-mAb-12584
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLRX5 C14orf87
Immunogen	Synthesized peptide derived from human glrx5
Specificity	This antibody detects endogenous levels of glrx5 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.76% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRX5;Glutaredoxin-related protein 5;mitochondrial;GLRX5;GLRX 5;FLB4739; GLRX5 C14orf87;mitochondrial (Monothiol glutaredoxin-5); mitochondrial (Monothiol glutaredoxin-5)
Observed Band	10-15 kDa
Calculated Molecular Weight	157 aa, 17 kDa
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images