



TMEM39A mouse mAb

Catalog No	YP-mAb-12585
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM39A
Immunogen	Synthesized peptide derived from human TMEM39A
Specificity	This antibody detects endogenous levels of TMEM39A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.77% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SUSR2;TMEM39A;transmembrane protein 39A
Observed Band	48-56 kDa
Calculated Molecular Weight	488 aa, 56 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Up-regulated in brain tumor glioblastoma multiforme cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images