



COMMD6 mouse mAb

Catalog No	YP-mAb-12590
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COMMD6 MSTP076
Immunogen	Synthesized peptide derived from human COMMD6 AA range: 1-50
Specificity	This antibody detects endogenous levels of COMMD6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.82% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COMMD6;COMMD6 MSTP076;COMM domain-containing protein 6
Observed Band	
Calculated Molecular Weight	9KD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Ubiquitous. Expressed in brain, heart, skeletal muscle, lung, pancreas, liver, kidney, small intestine and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images