



NeuN Neuronal Marker mouse mAb

Catalog No	YP-mAb-12591
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBFOX3
Immunogen	Synthesized peptide derived from human NeuN Neuronal Marker
Specificity	This antibody detects endogenous levels of NeuN Neuronal Marker at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.83% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fox-1 homolog C;FOX3;HRNBP3;NeuN antigen;Neuronal nuclei antigen;NeuN; RBFOX3;NeuN Neuronal Marker;RNA binding protein fox-1 homolog 3 (Fox-1 homolog C)
Observed Band	46-52 kDa
Calculated Molecular Weight	46 48kD
Cell Pathway	Nucleus . Cytoplasm . Largely restricted to neuronal nuclei. However, significant cytoplasmic localization in neurons from brains from HIV-infected individuals with cognitive impairment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images