



LDLRAD3 Monoclonal Antibody

Catalog No	YP-mAb-13392
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LDLRAD3
Immunogen	The antiserum was produced against synthesized peptide derived from human LDLRAD3. AA range:37-86
Specificity	LDLRAD3 Monoclonal Antibody detects endogenous levels of LDLRAD3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LDLRAD3;Low-density lipoprotein receptor class A domain-containing protein 3
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed at high levels in brain, lung, skeletal muscle, and pancreas. Expressed at moderate levels in heart, placenta, and kidney but not detected in the liver.
Function	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,
Background	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images