



CD92 Monoclonal Antibody

Catalog No	YP-mAb-13917
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC44A1
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC44A1. AA range:391-440
Specificity	CD92 Monoclonal Antibody detects endogenous levels of CD92 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD92;CDw92;Choline transporter-like protein 1;CHTL1;CTL1;SLC44A1;Solute carrier family 44 member 1;CD antigen CD92
Observed Band	70-73 kDa
Calculated Molecular Weight	73 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Mitochondrion outer membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in various cells of the hematopoietic system.
Function	function:Probable choline transporter. May be involved in membrane synthesis and myelin production.,similarity:Belongs to the CTL (choline transporter-like) family.,tissue specificity:Expressed in various cells of the hematopoietic system.,
Background	function:Probable choline transporter. May be involved in membrane synthesis and myelin production.,similarity:Belongs to the CTL (choline transporter-like) family.,tissue specificity:Expressed in various cells of the hematopoietic system.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images