



PI 3-kinase p101 Monoclonal Antibody

Catalog No	YP-mAb-14912
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PIK3R5
Immunogen	The antiserum was produced against synthesized peptide derived from human PIK3R5. AA range:695-744
Specificity	PI 3-kinase p101 Monoclonal Antibody detects endogenous levels of PI 3-kinase p101 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5-bisphosphate 3-kinase regulatory subunit;F7338I15Rik;FOAP 2;FOAP2; OTTHUMP135339;p11;P11 PI3K;p11-PI3K;Phosphatidylinositol 4;5 bisphosphate 3 kinase regulatory subunit;Phosphoinositide 3 kinase;regulatory subunit 5;regulatory subunit;polypeptide p11;Phosphoinositide 3-kinase regulatory subunit 5;PI3 kinase p11 subunit;PI3 kinase regulatory subunit 5;PI3-kinase p11 subunit;PI3-kinase regulatory subunit 5;PI3R5_HUMAN;PIK3 R5; PIK3R5;Protein FOAP 2;Protein FOAP-2;PtdIns 3 kinase p11;PtdIns 3 kinase regulatory subunit;PtdIns-3-kinase p11;PtdIns-3-kinase regulatory subunit;PI 3-kinase p101;PI3-kinase p101 subunit;PtdIns-
Observed Band	100kD
Calculated Molecular Weight	100kD
Cell Pathway	Nucleus . Cytoplasm . Cell membrane ; Peripheral membrane protein . Predominantly localized in the nucleus in absence of PIK3CG/p120. Colocalizes with PIK3CG/p120 in the cytoplasm. Translocated to the plasma membrane in a beta-gamma G protein-dependent manner. .
Tissue Specificity	Ubiquitously expressed with high expression in fetal brain compared to adult brain. Abundant expression is observed in cerebellum, cerebral cortex, cerebral meninges, and vermis cerebelli.



Function	domain:The heterodimerization region allows the binding to the catalytic subunit., enzyme regulation:Greatly activated by G gamma proteins.,function:Regulatory subunit of the PI3K gamma complex.,subunit:Heterodimer of a catalytic subunit (PIK3CG/p120) and a regulatory (PIK3R5a/p101) subunit. Interacts with G beta gamma proteins.,tissue specificity:Highly expressed in leukocytes, followed by spleen, lymph node, thymus and bone marrow.,
Background	Phosphatidylinositol 3-kinases (PI3Ks) phosphorylate the inositol ring of phosphatidylinositol at the 3-prime position, and play important roles in cell growth, proliferation, differentiation, motility, survival and intracellular trafficking. The PI3Ks are divided into three classes: I, II and III, and only the class I PI3Ks are involved in oncogenesis. This gene encodes the 101 kD regulatory subunit of the class I PI3K gamma complex, which is a dimeric enzyme, consisting of a 110 kD catalytic subunit gamma and a regulatory subunit of either 55, 87 or 101 kD. This protein recruits the catalytic subunit from the cytosol to the plasma membrane through high-affinity interaction with G-beta-gamma proteins. Multiple alternatively spliced transcript variants encoding two distinct isoforms have been found. [provided by RefSeq, Oct 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images