



EPG5 Monoclonal Antibody

Catalog No	YP-mAb-15054
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EPG5 KIAA1632
Immunogen	Synthetic peptide from human protein at AA range: 141-190
Specificity	The antibody detects endogenous EPG5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEEW1;KIAA1632;UPF0493 protein KIAA1632;EPG5;EPG5 KIAA1632;Ectopic P granules protein 5 homolog
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, perinuclear region . Lysosome .
Tissue Specificity	Brain,Cervix,Lymph node,Placenta,
Function	sequence caution:Translated as Trp.,similarity:Belongs to the UPF0493 family.,
Background	This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions. Mutations in this gene cause Vici syndrome. [provided by RefSeq, Aug 2015],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images