



Thyroperoxidase mouse mAb

Catalog No	YP-mAb-15720
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPO
Immunogen	Synthesized peptide derived from human Thyroperoxidase AA range: 651-700
Specificity	This antibody detects endogenous levels of Human Thyroperoxidase
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MSA;TPX;TDH2A;Thyroid Peroxidase;TPO;Thyroperoxidase;Thyroid peroxidase (TPO;EC 1.11.1.8)
Observed Band	103kD
Calculated Molecular Weight	103kD
Cell Pathway	Membrane; Single-pass type I membrane protein.; [Isoform 3]: Cell surface.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images