



PDCD1 mouse mAb

Catalog No	YP-mAb-15725
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDCD1 PD1
Immunogen	Synthesized peptide derived from human PDCD1
Specificity	This antibody detects endogenous levels of Human PDCD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDCD1 PD1;PDCD1;Programmed cell death protein 1 (Protein PD-1;hPD-1;CD antigen CD279)
Observed Band	31680DA
Calculated Molecular Weight	31680
Cell Pathway	Cell membrane ; Single-pass type I membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images