



ITIH4 (80k, Cleaved-Arg689) mouse mAb

Catalog No	YP-mAb-15818
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ITIH4 IHRP ITIHL1 PK120 PRO1851
Immunogen	Synthesized peptide derived from human ITIH4 (80k, Cleaved-Arg689)
Specificity	This antibody detects endogenous levels of Human ITIH4 (80k, Cleaved-Arg689)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	35 kDa inter-alpha-trypsin inhibitor heavy chain H4;70 kDa inter-alpha-trypsin inhibitor heavy chain H4;IHRP;Inter-alpha-inhibitor heavy chain 4;Inter-alpha-trypsin inhibitor family heavy chain-related protein;ITIH4;ITIH4 IHRP ITIHL1 PK120 PRO1851;ITIH4 (80k;Cleaved-Arg689);Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4;ITI-HC4;Plasma kallikrein sensitive glycoprotein 120;Gp120;PK-120) [Cleaved into: 70 kDa inter-alpha-trypsin inhibitor heavy chain H4;35 kDa inter-alpha-trypsin inhibitor heavy chain H4]; Cleaved-Arg689)
Observed Band	80 100kD
Calculated Molecular Weight	80 100kD
Cell Pathway	Secreted.
Tissue Specificity	Liver specific.
Function	#REF!
Background	#REF!
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images