



Kv1.5 Monoclonal Antibody

Catalog No	YP-mAb-16314
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCNA5
Immunogen	Synthetic Peptide of Kv1.5
Specificity	Kv1.5 protein(A257) detects endogenous levels of Kv1.5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCNA5;Potassium voltage-gated channel subfamily A member 5 (HPCN1) (Voltage-gated potassium channel HK2) (Voltage-gated potassium channel subunit Kv1.5);Potassium voltage-gated channel subfamily A member 5;HPCN1; Voltage-gated potassium channel HK2;Voltage-gated potassium channel subunit Kv1.5;YM3574
Observed Band	67
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Pancreatic islets and insulinoma.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images