



Cadherin-24 Monoclonal Antibody

Catalog No	YP-mAb-16933
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CDH24
Immunogen	The antiserum was produced against synthesized peptide derived from human CDH24. AA range:81-130
Specificity	Cadherin-24 Monoclonal Antibody detects endogenous levels of Cadherin-24 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cadherin-24;Cadherin 24;cadherin like 24;CDH11L;UNQ2834/PRO34009;CDH24
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Blood,Colon,HeLa cells,Placenta,Salivary gland,Testis,Thyroid,
Function	function:Cadherins are calcium dependent cell adhesion proteins., function:Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion.,similarity:Contains 5 cadherin domains.,subunit:Associates with alpha-, beta- and delta-catenins.,
Background	function:Cadherins are calcium dependent cell adhesion proteins., function:Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting



cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion.,similarity:Contains 5 cadherin domains.,subunit:Associates with alpha-, beta- and delta-catenins.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images