



NUPR1 mouse mAb

Catalog No	YP-mAb-17235
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NUPR1 COM1
Immunogen	Synthesized peptide derived from human C-terminal NUPR1
Specificity	This antibody detects endogenous levels of NUPR1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COM1;Nuclear protein 1;NUPR 1;NUPR1;NUPR1 COM1;Nuclear protein 1 (Candidate of metastasis 1) (Protein p8)
Observed Band	9kD
Calculated Molecular Weight	8 kDa
Cell Pathway	Nucleus . Cytoplasm . Cytoplasm, perinuclear region .
Tissue Specificity	Widely expressed, with high levels in liver, pancreas, prostate, ovary, colon, thyroid, spinal cord, trachea and adrenal gland, moderate levels in heart, placenta, lung, skeletal muscle, kidney, testis, small intestine, stomach and lymph node, and low levels in brain, spleen, thymus and bone marrow. Not detected in peripheral blood leukocytes.
Function	Transcription regulator that converts stress signals into a program of gene expression that empowers cells with resistance to the stress induced by a change in their microenvironment. Thereby participates in regulation of many process namely cell-cycle, apoptosis, autophagy and DNA repair responses . Controls cell cycle progression and protects cells from genotoxic stress induced by doxorubicin through the complex formation with TP53 and EP300 that binds CDKN1A promoter leading to transcriptional induction of CDKN1A . Protects pancreatic cancer cells from stress-induced cell death by binding the RELB



promoter and activating its transcription, leading to IER3 transactivation .
Negatively regulates apoptosis through interaction with PTMA . Inhibits
autophagy-induced apoptosis in cardiac cells through FOXO3 interaction,
inducing cytoplasmic translocation of FOXO3 thereby preventing the FOX

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images